

Cytology Test Listing

Test Code	Test	Specimen Requirements
104 or 1104	Pap Smear	A properly collected, fixed and labeled specimen - See Specimen Collection Techniques for conventional pap smears for details. A cytology requisition that is completely and properly filled out – See Cytology Requisition for details Place specimen and requisition in transport bag and place in area designated for courier pickup. – See Specimen Transport of Cytology Specimens for details.
108 or 1108	Thin Layer Pap Smear	A properly collected, fixed and labeled specimen - See Specimen Collection Techniques for thin layer pap smears for details. A cytology requisition that is completely and properly filled out – See Cytology Requisition for details Place specimen and requisition in transport bag and place in area designated for courier pickup. – See Specimen Transport of Cytology Specimens for details.
137	Sputum	A properly collected, fixed and labeled specimen. - See Specimen Collection Techniques for sputum for details. A cytology requisition that is completely and properly filled out – See Cytology Requisition for details Place specimen and requisition in transport bag and place in area designated for courier pickup. – See Specimen Transport of Cytology Specimens for details
172	Bronchial Washings	A properly collected, fixed and labeled specimen - See Specimen Collection Techniques for washings for details. A cytology requisition that is completely and properly filled out – See Cytology Requisition for details Place specimen and requisition in transport bag and place in area designated for courier pickup. – See Specimen Transport of Cytology Specimens for details
175	Bronchial Brushing	A properly collected, fixed and labeled specimen - See Specimen Collection Techniques for brushing specimens for details. A cytology requisition that is completely and properly filled out – See Cytology Requisition for details Place specimen and requisition in transport bag and place in area designated for courier pickup. – See Specimen Transport of Cytology Specimens for details
290	Fine Needle Aspirations	A properly collected, fixed and labeled specimen - See Specimen Collection Techniques for FNA (solid and cystic) specimens for details. A cytology requisition that is completely and properly filled out, indicate exact specimen site. Example: FNA of thyroid not FNA of neck. See Cytology Requisition for details. Place specimen and requisition in transport bag and place in area designated for courier pickup. – See Specimen Transport of Cytology Specimens for details
110	Body Fluids	Includes large/small volume fluids, non-bronchial brushings/washings, urines, herpes smears. A properly collected, fixed and labeled specimen - See Specimen Collection Techniques for body fluids for details. A cytology requisition that is completely and properly filled out – See Cytology Requisition for details Place specimen and requisition in transport bag and place in area designated for courier pickup. – See Specimen Transport of Cytology Specimens for details
125	HPV Typing	Additional test that may be ordered on thin layer pap smears. Indicate in the available area on the requisition.
127	Chlamydia	Additional test that may be ordered on thin layer pap smears. Indicate in the “other” area of the requisition.
129	Neisseria Gonorrhoeae	Additional test that may be ordered on thin layer pap smears. Indicate in the “other” area of the requisition.